

Package ‘sRm’

November 14, 2017

Type Package

Title Minimal parser for Selective Reaction Monitoring (SRM) mass spectrometry (MS) data

Version 0.1.4

Author Tom Wilson <tpw2@aber.ac.uk>

Maintainer Tom Wilson <tpw2@aber.ac.uk>

Description Functions for parsing SRM-MS data from .mzML files.

License GPL (>= 3)

Depends methods, xml2 (>= 1.0), ggplot2, grid, gridExtra, dplyr

Imports mzR (>= 2.1), ggrepel, tibble, purrr

URL <https://github.com/wilsontom/sRm>

BugReports <https://github.com/wilsontom/sRm/issues>

LazyData TRUE

RoxygenNote 6.0.1

Suggests testthat

NeedsCompilation no

R topics documented:

bpi	2
combineTransitions	2
export_to_skyline	3
get	4
get_meta	4
get_scan_header	5
meta	6
openSRMfile	6
plotAll	7
plotMulti	7
plotSRM	8

show,SRM-method	8
SRM-class	9
tic	9
transitions	10

Index	11
--------------	-----------

bpi	<i>bpi</i>
-----	------------

Description

bpi

Usage

```
bpi(object)

## S4 method for signature 'SRM'
bpi(object)
```

Arguments

object a SRM object

Value

a numeric value for the whole sample base peak intensity (bpi)

Author(s)

Tom Wilson <tpw2@aber.ac.uk>

combineTransitions	<i>Combine Transitions</i>
--------------------	----------------------------

Description

Combine Transitions

Each parent m/z will usually have multiple separate Q3 product ions. SRM transitions are stored separately in the mzML file. Each index is a transition from a single Q1 precursor ion to one of the specified Q3 ions. This method combines all Q3 product ion chromatograms into one, to produce a total ion chromatogram based on the parent m/z and all Q3 product ions assigned to that parent m/z

Where a dynamic MRM is used; then the start retention time (Rt) and end Rt is also taken into account when transitions are combined.

For example the following transitions from a standard SRM-MS experiment would be combined into a single chromatogram based on the parent m/z value;

- SRM SIC 153.01,65.271
- SRM SIC 153.01,67.232
- SRM SIC 153.01,109.094

Usage

```
combineTransitions(object)
```

```
## S4 method for signature 'SRM'  
combineTransitions(object)
```

Arguments

object a SRM object

Value

a transition class of SRM total ion chromatograms

Author(s)

Tom Wilson <tpw2@aber.ac.uk>

export_to_skyline	<i>export_to_skyline</i>
-------------------	--------------------------

Description

export_to_skyline

Usage

```
export_to_skyline(object, transitions = NULL)
```

```
## S4 method for signature 'SRM'  
export_to_skyline(object, transitions = NULL)
```

Arguments

object	a SRM object
transitions	an optional data.frame containing the following columns for selected (or all) transitions.

Value

a tibble of transitions in a format which can be easily used with Skyline

Author(s)

Tom Wilson <tpw2@aber.ac.uk>

get	get
-----	-----

Description

get

Usage

```
get(object, Q1, Q3, method)
```

```
## S4 method for signature 'SRM'
```

```
get(object, Q1, Q3, method)
```

Arguments

object	a SRM object
Q1	a numeric value for parent m/z
Q3	a numeric value for product m/z
method	either 'tic' or 'bpi' to return total ion count or base peak intensity respectively

Value

a tibble of parent and product mass, and either tic or bpi

get_meta	Extract meta information from .mzML file
----------	--

Description

Extract a minimal amount meta data from a .mzML file

Usage

```
get_meta(x)
```

Arguments

x	a valid .mzML file
---	--------------------

Value

a tibble containing;

- mzML Schema
- Acquisition Date
- Acquisition Time
- Instrument Model
- File ID

get_scan_header	<i>Extract scan header information</i>
-----------------	--

Description

Extract key values (parent m/z, product m/z and polarity) from the chromatogram blocks of an .mzML file

Usage

```
get_scan_header(x)
```

Arguments

x a valid .mzML file

Value

a tibble containing

- header
- polarity
- Q1
- Q3
- tidy_head

meta	<i>Show meta data</i>
------	-----------------------

Description

Show meta data

Usage

```
meta(object)
```

```
## S4 method for signature 'SRM'  
meta(object)
```

Arguments

object a SRM object

Author(s)

Tom Wilson <tpw2@aber.ac.uk>

openSRMfile	<i>Open a SRM-MS file</i>
-------------	---------------------------

Description

Open and parse mzML file into a SRM object

Usage

```
openSRMfile(mzMLFile)
```

Arguments

mzMLFile a valid .mzML file

Value

a SRM object (see [SRM-class](#))

Author(s)

Tom Wilson <tpw2@aber.ac.uk>

plotAll	<i>Plot all</i>
---------	-----------------

Description

Plot the SRM-MS TIC and all individual Q3 traces for a selected parent scans'

Usage

```
plotAll(x, n)

## S3 method for class 'transition'
plotAll(x, n)
```

Arguments

x	a transition class
n	the number of the transition (index) to plot

Author(s)

Tom Wilson <tpw2@aber.ac.uk>

plotMulti	<i>plotMulti</i>
-----------	------------------

Description

plotMulti

Usage

```
plotMulti(object, idn, addLabels = FALSE, labels = NULL)

## S4 method for signature 'SRM'
plotMulti(object, idn, addLabels = TRUE, labels = NULL)
```

Arguments

object	a SRM object
idn	the index number of a transition to plot
addLabels	logical; If TRUE then labels for each geom_line are added
labels	an optional character vector of label names. Default is NULL. If addLabels is TRUE and labels is NULL then object@index is used for label text.

Author(s)

Tom Wilson <tpw2@aber.ac.uk>

plotSRM

Plot SRM

Description

Plot SRM

Usage

```
plotSRM(object, idn)
```

```
## S4 method for signature 'SRM'
```

```
plotSRM(object, idn)
```

Arguments

object a SRM object

idn the index number of a transition to plot

Author(s)

Tom Wilson <tpw2@aber.ac.uk>

show,SRM-method

show-SRM

Description

show-SRM

Usage

```
## S4 method for signature 'SRM'
```

```
show(object)
```

Arguments

object a SRM object

Author(s)

Tom Wilson <tpw2@aber.ac.uk>

SRM-class

SRM

Description

A S4 class to store single reaction monitoring (SRM) mass spectrometry (MS) data

Slots

SHA1 character vector of the original SHA-1 checksum value

meta list containing file meta data

totIonCount data.frame of sample total ion count (TIC)

filter character vector of scan filters; taken directly from chromatogram idRefs

index character vector of readable scan index's; ie Q1: 153.01 --> Q3: 65.271 (-)

peaks list of peak data

header data.frame of summary information for each scan index;parent, product, polarity,
totIonCount, basePeakInt

tic

tic

Description

tic

Usage

```
tic(object)
```

```
## S4 method for signature 'SRM'
```

```
tic(object)
```

Arguments

object a SRM object

Value

a numeric value for the whole sample total ion count (tic)

Author(s)

Tom Wilson <tpw2@aber.ac.uk>

transitions

Show transitions

Description

Show transitions

Usage

```
transitions(object)
```

```
## S4 method for signature 'SRM'  
transitions(object)
```

Arguments

object a SRM object

Author(s)

Tom Wilson <tpw2@aber.ac.uk>

Index

bpi, [2](#)
bpi, SRM-method (bpi), [2](#)

combineTransitions, [2](#)
combineTransitions, SRM-method
(combineTransitions), [2](#)

export_to_skyline, [3](#)
export_to_skyline, SRM-method
(export_to_skyline), [3](#)

get, [4](#)
get, SRM-method (get), [4](#)
get_meta, [4](#)
get_scan_header, [5](#)

meta, [6](#)
meta, SRM-method (meta), [6](#)

openSRMfile, [6](#)

plotAll, [7](#)
plotMulti, [7](#)
plotMulti, SRM-method (plotMulti), [7](#)
plotSRM, [8](#)
plotSRM, SRM-method (plotSRM), [8](#)

show, SRM-method, [8](#)
SRM-class, [9](#)

tic, [9](#)
tic, SRM-method (tic), [9](#)
transitions, [10](#)
transitions, SRM-method (transitions), [10](#)